

MEDICAL CART



SPECIFICATION

Top & Bottom shelf : ABS plastic with SUS 304
Guard rail : Steel (Color : Grey)
Column : Aluminum (Color : Grey)
Drawer front : ABS plastic
Drawer with removable ABS basket



MEDICAL CART


Model	Size (L×D×H / mm)	Opstion
KS-MDC-D1	740×575×900	H:1000
KS-MDCS-D1	610×575×900	H:1000

■ OPTION 높이

MEDICAL CART

Model	Size (L×D×H / mm)	Opstion
KS-MDC-D1-P1	740×575×900	H:1000
KS-MDC-D1S-P1	610×575×900	H:1000

Applicable Components (Option)



MONITOR ARM

Model	Remark
KS-ARM-MSB	Height adjustable

■ VESA 75x75 / 100x100
■ OPTION



(New model)



BASKET

Model	Size (L×D×H / mm)
KS-AB-1	560×410×55
KS-AB-2	560×410×125
KS-AB-3	560×410×190

■ DIVIDER compartment control through
■ OPTION



WASTE BOX

Model	Size (L×D×H / mm)
KS-AC-WB01	185×175×330 (10L)

■ OPTION



IV POLE

Model	Size (L×D×H / mm)
KS-AC-IV05	φ22

■ H=450
■ 5 HOOKS
■ OPTION



WASTE HOLDER

Model	Size (L×D×H / mm)
KS-AC-WH02	210×230×155 (10L)
KS-AC-WH03	280×300×155 (20L)

■ Application of new model
■ OPTION

NEEDLE HOLDER

Model	Size (L×D×H / mm)
KS-AC-NH02	155×120×110 (2L)
KS-AC-NH03	200×150×110 (5L)

내측사이즈



MEDICAL CART

Model	Size (L×D×H / mm)	Opstion
KS-MDC	740×575×900	H:1000
KS-MDCS	610×575×900	H:1000



MEDICAL CART

Model	Size (L×D×H / mm)	Opstion
KS-MDC-YD1	740×575×900	H:1000
KS-MDCS-YD1	610×575×900	H:1000



MEDICAL CART

Model	Size (L×D×H / mm)	Opstion
KS-MDC-P1	740×575×900	H:1000
KS-MDCS-P1	610×575×900	H:1000



MEDICAL CART

Model	Size (L×D×H / mm)	Opstion
KS-MDC-D2	740×575×900	H:1000
KS-MDCS-D2	610×575×900	H:1000



MEDICAL CART

Model	Size (L×D×H / mm)	Opstion
KS-MDCS-D2W2	740×575×900	H:1000



MEDICAL CART

Model	Size (L×D×H / mm)	Opstion
KS-MDC-D3	740×575×1000	H:1000
KS-MDCS-D3	610×575×1000	H:1000



MEDICAL CART

Model	Size (L×D×H / mm)	Opstion
KS-MDC-KD1-P1	740×575×900	H:1000
KS-MDCS-KD1-P1	610×575×900	H:1000

■ Keyboard



MEDICAL CART

Model	Size (L×D×H / mm)	Opstion
KS-MDC-KD2	740×575×900	H:1000
KS-MDCS-KD2	610×575×900	H:1000

■ Keyboard



IV POLE (H 2단행거)

Model	Size (L×D×H / mm)	Opstion
KS-AC-IV06	φ22	H:1000

EMERGENCY CART




Model	KS-EMC5D-F
-------	------------



ALUMINUM SIDE & BACK PANEL
Light weight of ALUMINUM makes faster moving speed than cart made by steel.

SPECIFICATION
Top shelf : ABS plastic with SUS 304
Side & Back panel : Aluminum
Guard rail : Steel (Color : Grey)
Column : Aluminum (Color : Grey)
Drawer front : ABS plastic
Drawer with removable ABS basket
Drawer rail : Auto return available upon request
Total key locking
Caster : 5 inch (Two with brake)
Height can be modified upon request



Model	Remark
KS-ARM-MSB	Height adjustable
■ VESA 75x75 / 100x100	
■ OPTION	



You can open the drawer and take out the basket.
You can move archived content.




Model	Size (L×D×H / mm)
KS-AC-OT01	For φ120

Oxygen tank, C.P.R Board can be fitted with. (Option)




Model	Size (L×D×H / mm)
KS-AC-CPR01	420×530×20

ACCESSORIES



Model	Size (L×D×H / mm)
KS-AC-DF01	L : 395X275 S : 330X235

OPTION



Model	Size (L×D×H / mm)
KS-AC-UC01	340X100X300(250)

OPTION



Model	Size (L×D×H / mm)
KS-AC-NH02	155×120×110 (2L)
KS-AC-NH03	200×150×110 (5L)

내측사이즈



Model	Size (L×D×H / mm)
-	350X340

SLIDE PLATE



Model	Size (L×D×H / mm)
KS-AC-SB01	480X180X418

6 PERSON
OPTION
KS-CA703



Model	Size (L×D×H / mm)
KS-AC-WB01	185×175×330

OPTION



Model	Size (L×D×H / mm)
KS-AB-1	560×410×55
KS-AB-2	560×410×125
KS-AB-3	560×410×190

DIVIDER Compartment adjustment
OPTION



Model	Size (L×D×H / mm)
KS-AC-IV05	φ22

H=450
5 HOOKS
OPTION



Model	Size (L×D×H / mm)
KS-AC-WH02	210×230×155 (10L)
KS-AC-WH03	280×300×155 (20L)

New model applied
OPTION

W740 SERIES



EMERGENCY CART (740)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D	740×575×985(1055)	Standard

- Lock (Key)



EMERGENCY CART (740)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-01	740×575×985(1055)	—

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)



EMERGENCY CART (740)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-02	740×575×985(1055)	—

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)



EMERGENCY CART (740)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-03	740×575×985(1055)	—

- Storage Box (6 person) 2ea
- Lock (Key)

W610 SERIES



EMERGENCY CART (610)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-S	610×575×985(1055)	Standard



EMERGENCY CART (610)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-01S	610×575×985(1055)	—

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder



EMERGENCY CART (610)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-02S	610×575×985(1055)	—

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder



EMERGENCY CART (610)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-03S	610×575×985(1055)	—

- Storage Box (6 person) 2ea

Operating Room
Appliances

Nursing
Appliances

Consulting Room
Appliances

Supply Room
Appliances

Hospital Bed

Carts

Monitor Carts

IV Pole Stand

Monitor Arm

Wall Mount


Pole Mount

Accessory

MEDICATION CART (BIN TYPE)



ACCESSORIES



DEFIBRILLATOR

Model	Size (L×D×H / mm)
KS-AC-DF01	L : 395X275 S : 330X235

OPTION




SLIDE PLATE

Model	Size (L×D×H / mm)
KS-AC-SP1	350X340



NEEDLE HOLDER

Model	Size (L×D×H / mm)
KS-AC-SB02	480×180×418



UTILITY CONTAINER

Model	Size (L×D×H / mm)
KS-AC-UC01	340X100X300(250)

OPTION



STORAGE BOX

Model	Size (L×D×H / mm)
KS-AC-SB01	480X180X418

6 PERSON
OPTION
KS-CA703



IV POLE

Model	Size (L×D×H / mm)
KS-AC-IV05	φ22

H=450
5 HOOKS
OPTION

Model	KS-MCC30-F
-------	------------

AVAILABLE BIN TYPE

- 16 bins type (Front side 8 bins + Back side 8 bins)
- 20 bins type (Front side 10 bins + Back side 10 bins)
- 24 bins type (Front side 12 bins + Back side 12 bins)
- 30 bins type (Front side 15 bins + Back side 15 bins)
- 32 bins type (Front side 16 bins + Back side 16 bins)
- 40 bins type (Front side 20 bins + Back side 20 bins)

ALUMINUM SIDE & BACK PANEL


Light weight of ALUMINUM makes faster moving speed than cart made by steel.

SPECIFICATION

SIZE (L×D×H / mm) : 740 X 575 X 1000(1070)
610 X 575 X 1000(1070)
Top shelf : ABS plastic with SUS 304
Side & Back panel : Aluminum
Guard rail : Steel (Color : Grey)
Column : Aluminum (Color : Grey)
Drawer : Steel with fixed one divider
Bin : ABS plastic
Total key locking
Caster : 5 inch (Two with brake)
Height can be modified upon request

MEDICATION BIN SIZE

KS-AC-MB01 (L×D×H / mm) : 108 X 248 X 86
ABS plastic
Including 2 removable dividers



You can adjust the space inside of the Bin Box by using a partition



TOTAL KEY LOCKING

Every Medication bins and Drawers at front & back side can be locked by one key. Use the key to fully lock the bin box to prepare for loss



NEEDLE HOLDER

Model	Size (L×D×H / mm)
KS-AC-NH02	155×120×110 (2L)
KS-AC-NH03	200×150×110 (5L)

내측사이즈



WASTE BOX

Model	Size (L×D×H / mm)
KS-AC-WB01	185×175×330

OPTION



WASTE HOLDER

Model	Size (L×D×H / mm)
KS-AC-WH02	210×230×155 (10L)
KS-AC-WH03	280×300×155 (20L)

New model applied
OPTION

W740 SERIES

Nursing
Appliances



MEDICATION CART (740,DRAWER) MEDICATION CART (740,DRAWER) MEDICATION CART (740,DRAWER)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC20	740×575×1000(1070)	20 Person	KS-MCC30	740×575×1000(1070)	30 Person	KS-MCC40	740×575×1000(1070)	40 Person

- Standard
■ Lock (Key)
- Standard
■ Lock (Key)
- Standard
■ Lock (Key)



MEDICATION CART (740,DRAWER) MEDICATION CART (740,DRAWER) MEDICATION CART (740,DRAWER)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC20-01	740×575×1000(1070)	20 Person	KS-MCC30-01	740×575×1000(1070)	30 Person	KS-MCC40-01	740×575×1000(1070)	40 Person

- Utility Container
■ Waste Holder
■ Waste Box
■ Needle Holder
■ Lock (Key)
- Utility Container
■ Waste Holder
■ Waste Box
■ Needle Holder
■ Lock (Key)
- Utility Container
■ Waste Holder
■ Waste Box
■ Needle Holder
■ Lock (Key)

W740 SERIES



MEDICATION CART (740,DRAWER) MEDICATION CART (740,DRAWER) MEDICATION CART (740,DRAWER)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC20-02	740×575×1000(1070)	20 Person	KS-MCC30-02	740×575×1000(1070)	30 Person	KS-MCC40-02	740×575×1000(1070)	40 Person

- Storage Box (6 person)
■ Waste Holder
■ Waste Box
■ Needle Holder
■ Lock (Key)
- Storage Box (6 person)
■ Waste Holder
■ Waste Box
■ Needle Holder
■ Lock (Key)
- Storage Box (6 person)
■ Waste Holder
■ Waste Box
■ Needle Holder
■ Lock (Key)



MEDICATION CART (740,DRAWER) MEDICATION CART (740,DRAWER) MEDICATION CART (740,DRAWER)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC20-03	740×575×1000(1070)	20 Person	KS-MCC30-03	740×575×1000(1070)	30 Person	KS-MCC40-03	740×575×1000(1070)	40 Person

- Storage Box (6 person) 2ea
■ Lock (Key)
- Storage Box (6 person) 2ea
■ Lock (Key)
- Storage Box (6 person) 2ea
■ Lock (Key)

Operating Room
Appliances
Nursing
Appliances
Consulting Room
Appliances
Supply Room
Appliances
Hospital Bed
Carts
Monitor Carts
IV Pole Stand
Monitor Arm
Wall Mount
Pole Mount
Accessory

W610 SERIES

Nursing
Appliances



MEDICATION CART (610,DRAWER)			MEDICATION CART (610,DRAWER)			MEDICATION CART (610,DRAWER)		
Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC16-01	610×575×1000(1070)	16 Person	KS-MCC24-01	610×575×1000(1070)	24 Person	KS-MCC32-01	610×575×1000(1070)	32 Person
■ Standard ■ Lock (Key)			■ Standard ■ Lock (Key)			■ Standard ■ Lock (Key)		



MEDICATION CART (610,DRAWER)			MEDICATION CART (610,DRAWER)			MEDICATION CART (610,DRAWER)		
Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC16-02	610×575×1000(1070)	16 Person	KS-MCC24-02	610×575×1000(1070)	24 Person	KS-MCC32-02	610×575×1000(1070)	32 Person
■ Utility Container ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key)			■ Utility Container ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key)			■ Utility Container ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key)		

W610 SERIES



MEDICATION CART (610,DRAWER)			MEDICATION CART (610,DRAWER)			MEDICATION CART (610,DRAWER)		
Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC16-02	610×575×1000(1070)	16 Person	KS-MCC24-02	610×575×1000(1070)	24 Person	KS-MCC32-02	610×575×1000(1070)	32 Person
■ Storage Box (6 person) ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key)			■ Storage Box (6 person) ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key)			■ Storage Box (6 person) ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key)		



MEDICATION CART (610,DRAWER)			MEDICATION CART (610,DRAWER)			MEDICATION CART (610,DRAWER)		
Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC16-03	610×575×1000(1070)	16 Person	KS-MCC24-03	610×575×1000(1070)	24 Person	KS-MCC32-03	610×575×1000(1070)	32 Person
■ Storage Box (6 person) 2ea ■ Lock (Key)			■ Storage Box (6 person) 2ea ■ Lock (Key)			■ Storage Box (6 person) 2ea ■ Lock (Key)		

Operating Room
Appliances

Nursing
Appliances

Consulting Room
Appliances

Supply Room
Appliances

Hospital Bed

Carts

Monitor Carts

IV Pole Stand

Monitor Arm

Wall Mount

Pole Mount

Accessory


MEDICATION CART (CASSETTE TYPE)



Model	KS-MCC10C-F
-------	-------------

AVAILABLE CASSETTE TYPE

- 8 bins / 12 bins / 16 bins type
- 10 bins / 15 bins / 20 bins type



MEDICATION BIN SIZE


KS-AC-MB01 (LxDxH / mm) : 104 X 300 X 78
ABS plastic
Including 2 removable dividers

ALUMINUM SIDE & BACK PANEL

Light weight of ALUMINUM makes faster moving speed than cart made by steel.


SPECIFICATION

SIZE (LxDxH / mm) : 740 X 575 X 1000(1070)
610 X 575 X 1000(1070)
Top shelf : ABS plastic with SUS 304
Side & Back panel : Aluminum
Guard rail : Steel (Color : Grey)
Column : Aluminum (Color : Grey)
Drawer front : ABS plastic
Drawer with removable ABS basket
Drawer rail : Auto return available upon request
Bin : ABS plastic
Cassette : ALUMINUM
Caster : 5 inch (Two with brake)
Height can be modified upon request



MEDICATION CASSETTE IS TOTAL LOCKING BY KEY AND REMOVABLE.

ACCESSORIES



NOTEBOOK PLATE

Model	Size (LxDxH / mm)
KS-AC-DF01	395X275

OPTION



SLIDE PLATE


Model	Size (LxDxH / mm)
—	350X340



BASKET

Model	Size (LxDxH / mm)
KS-AB-1	560 X 410 X 55
KS-AB-2	560 X 410 X 125
KS-AB-3	560 X 410 X 190


DIVIDER compartment control through
OPTION



UTILITY CONTAINER

Model	Size (LxDxH / mm)
KS-AC-UC01	340X100X300(250)

OPTION



STORAGE BOX

Model	Size (LxDxH / mm)
KS-AC-SB01	480X180X418

6 PERSON
KS-CA703



IV POLE

Model	Size (LxDxH / mm)
KS-AC-IV05	φ22

H=450
5 HOOKS
OPTION



NEEDLE HOLDER

Model	Size (LxDxH / mm)
KS-AC-NH02	155 X 120 X 110 (2L)
KS-AC-NH03	200 X 150 X 110 (5L)

내측사이즈



WASTE BOX

Model	Size (LxDxH / mm)
KS-AC-WB01	185 X 175 X 330

OPTION



WASTE HOLDER

Model	Size (LxDxH / mm)
KS-AC-WH02	210 X 230 X 155 (10L)
KS-AC-WH03	280 X 300 X 155 (20L)

New model applied
OPTION

W740 SERIES

Nursing
Appliances



MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC10C	740×575×1000(1070)	10 Person	KS-MCC15C	740×575×1000(1070)	15 Person	KS-MCC20C	740×575×1000(1070)	20 Person
■ Standard ■ Lock (Key) ■ Cassette Type			■ Standard ■ Lock (Key) ■ Cassette Type			■ Standard ■ Lock (Key) ■ Cassette Type		



MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC10C-01	740×575×1000(1070)	10 Person	KS-MCC15C-01	740×575×1000(1070)	15 Person	KS-MCC20C-01	740×575×1000(1070)	20 Person
■ Utility Container ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key) ■ Cassette Type			■ Utility Container ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key) ■ Cassette Type			■ Utility Container ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key) ■ Cassette Type		

W740 SERIES



MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC10C-02	740×575×1000(1070)	10 Person	KS-MCC15C-02	740×575×1000(1070)	15 Person	KS-MCC20C-02	740×575×1000(1070)	20 Person
■ Storage Box (6 person) ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key) ■ Cassette Type			■ Storage Box (6 person) ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key) ■ Cassette Type			■ Storage Box (6 person) ■ Waste Holder ■ Waste Box ■ Needle Holder ■ Lock (Key) ■ Cassette Type		



MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC10C-03	740×575×1000(1070)	10 Person	KS-MCC15C-03	740×575×1000(1070)	15 Person	KS-MCC20C-03	740×575×1000(1070)	20 Person
■ Storage Box (6 person) 2ea ■ Lock (Key)			■ Storage Box (6 person) 2ea ■ Lock (Key)			■ Storage Box (6 person) 2ea ■ Lock (Key)		

Operating Room
Appliances
Nursing
Appliances
Consulting Room
Appliances
Supply Room
Appliances
Hospital Bed
Carts
Monitor Carts
IV Pole Stand
Monitor Arm
Wall Mount
Pole Mount
Accessory

W610 SERIES

Nursing
Appliances



MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC8C	610×575×1000(1070)	8 Person	KS-MCC12C	610×575×1000(1070)	12 Person	KS-MCC16C	610×575×1000(1070)	16 Person

- Standard
- Lock (Key)
- Cassette Type



MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC8C-01	610×575×1000(1070)	—	KS-MCC12C-01	610×575×1000(1070)	—	KS-MCC16C-01	610×575×1000(1070)	—

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type

W610 SERIES



MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC8C-02	610×575×1000(1070)	8 Person	KS-MCC12C-02	610×575×1000(1070)	12 Person	KS-MCC16C-02	610×575×1000(1070)	16 Person

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type



MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE)

Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark	Model	Size (L×D×H / mm)	Remark
KS-MCC8C-03	610×575×1000(1070)	8 Person	KS-MCC12C-03	610×575×1000(1070)	12 Person	KS-MCC16C-03	610×575×1000(1070)	16 Person

- Storage Box (6 person) 2ea
- Lock (Key)
- Cassette Type

Operating Room
Appliances

Nursing
Appliances

Consulting Room
Appliances

Supply Room
Appliances

Hospital Bed

Carts

Monitor Carts

IV Pole Stand

Monitor Arm

Wall Mount

Pole Mount

Accessory



MEDICATION CART

Model	Size (L×D×H / mm)	Remark
KS-SMCC25C-D1	740×575×1000	—



MEDICATION CART

Model	Size (L×D×H / mm)	Remark
KS-MC2	630×520×950	SUS(40인용)

■ 20/30 person Custom made is available.
■ CA502-20EA



CASSETTE DELIVERY CART

Model	Size (L×D×H / mm)	Remark
KS-CDLC-4B	740×740×1400	4 BOX



front and back



MEDICATION CART

Model	Size (L×D×H / mm)	Remark
KS-CDLC-4BL	740×740×1000	—



MEDICATION CART

Model	Size (L×D×H / mm)	Remark
KS-CDLC-4BS	740×740×1200	—



CASSETTE

Model	Size (L×D×H / mm)	Remark
KS-CB-10	563×330×250	10 Person

■ Aluminium, A.B.S Box



CASSETTE

Model	Size (L×D×H / mm)	Remark
KS-CB-15	563×330×345	15 Person

■ Aluminium, A.B.S Box



CASSETTE

Model	Size (L×D×H / mm)	Remark
KS-CB-20	563×330×440	20 Person

■ Aluminium, A.B.S Box



CASSETTE

Model	Size (L×D×H / mm)	Remark
KS-CB-8	463×330×250	8 Person

■ Aluminium, A.B.S Box



CASSETTE

Model	Size (L×D×H / mm)	Remark
KS-CB-12	463×330×345	12 Person

■ Aluminium, A.B.S Box



CASSETTE

Model	Size (L×D×H / mm)	Remark
KS-CB-16	463×330×440	16 Person

■ Aluminium, A.B.S Box



BIN BOX (BIN TYPE)

Model	Size (L×D×H / mm)	Remark
KS-AC-MB01	108×248×86	A.B.S Plastic

■ Space can be managed by adjusting the space.




BIN BOX (CASSETTE TYPE)

Model	Size (L×D×H / mm)	Remark
KS-AC-MB02	104×300×78	A.B.S Plastic

■ Space can be managed by adjusting the space.

ANESTHESIA CART



ALUMINUM SIDE & BACK PANEL

Light weight of ALUMINUM makes faster moving speed than cart made by steel.

SPECIFICATION

SIZE (LxDxH / mm) : 740 X 575 X 1070(1670)
610 X 575 X 1070(1670)

Top shelf : ABS plastic with SUS 304

Side & Back panel : Aluminum

Guard rail : Steel (Color : Grey)

Column : Aluminum (Color : Grey)

Drawer front : ABS plastic

Drawer with removable ABS basket

Drawer rail : Auto return available upon request

Multi Bin : ABS plastic (4, 5, 6 BIN available)

Container shelf : Steel

Total key locking

Caster : 5 inch (Two with brake)

Height can be modified upon request

ANESTHESIA CART

Model	Size (LxDxH / mm)	Remark
KS-ATC-01-01	740 X 575 X 1070(1670)	—



ANESTHESIA CART

Model	Size (LxDxH / mm)	Remark
KS-ATC-01-02	740 X 575 X 1070(1670)	—



ANESTHESIA CART

Model	Size (LxDxH / mm)	Remark
KS-ATC-01-03	740 X 575 X 1070(1670)	10 Hooks



ICU CART

Model	Size (LxDxH / mm)	Remark
KS-ICUC-01	740 X 575 X 1000(1070)	—



SUPER PLUS CART

Model	Size (LxDxH / mm)	Remark
KS-TRTC-5D	755X530X1060	Door Type



UTILITY CART

Model	Size (LxDxH / mm)	Remark
KS-ANC-02	660X500X1000	—



ANESTHESIA CART

Model	Size (LxDxH / mm)	Remark
KS-ANC-01	660X500X(1000+600)	—



SUPER PLUS CART

Model	Size (LxDxH / mm)	Remark
KS-TRTC-5D-01	755X530X1060	Open Type



UTILITY CART

Model	Size (LxDxH / mm)	Remark
KS-ANC-03	660X500X1060	—



WHEELCHAIR

Model	Size (L×D×H / mm)	Remark
KS-KCA1004	745X1150X875	Steel colossus



WHEELCHAIR

Model	Size (L×D×H / mm)	Remark
KS-KCA100	700X1040X875	Steele general



UTILITY CART (ABS DRAWER) KEY BAR

Model	Size (L×D×H / mm)	Remark
KS-ANC-04	660×500×1060	—



UTILITY CART

Model	Size (L×D×H / mm)	Remark
KS-ANC-05	660×500×1060	—



CO2 SUPPORT

Model	Size (L×D×H / mm)	Remark
KS-KCA100-02		



STAND HOLDER (RED CIRCLE)

Model	Size (L×D×H / mm)	Remark
KS-KCA100-PM		



UTILITY CART(DRAWER TYPE)

Model	Size (L×D×H / mm)	Remark
KS-UC3S-1D	665X500X950	



UTILITY CART (DRAWER TYPE)

Model	Size (L×D×H / mm)	Remark
KS-UC2S-2S	665X500X950	2 DRAWERS




WALKER

Model	Size (L×D×H / mm)	Remark
KS-WK01	750×450×950	STANDARD
KS-WK02	630X500X950	TWO BUTTON



HAIR SHAMPOOING BASKET

Model	Size (L×D×H / mm)	Remark
KS-HS-PA	450×440×85	Includes water cooker



UTILITY CART (DRAWER TYPE)

Model	Size (L×D×H / mm)	Remark
KS-UC25-BP	650X350X950	



UTILITY CART (A.B.S DRAWER)

Model	Size (L×D×H / mm)
KS-AUC2S-1DS	565 × 500 × 950
KS-AUC2S-1DL	665 × 500 × 950

Stainless steel 304 construction
3 side guard rail on each shelf
Two flat shelves
Drawer front : ABS plastic
ABS basket with divider each drawer
Round rubber bumper
4 inch casters (Two with breaks)

UTILITY CART (A.B.S DRAWER)

Model	Size (L×D×H / mm)
KS-AUC3S-1DS	565 × 500 × 950
KS-AUC3S-1DL	665 × 500 × 950

Stainless steel 304 construction
3 side guard rail on each shelf
Three flat shelves
Drawer front : ABS plastic
ABS basket with divider each drawer
Round rubber bumper
4 inch casters (Two with breaks)

UTILITY CART (A.B.S DRAWER)

Model	Size (L×D×H / mm)
KS-AUC2S-2DS	565 × 500 × 950
KS-AUC2S-2DL	665 × 500 × 950

Stainless steel 304 construction
3 side guard rail on each shelf
Two flat shelves
Drawer front : ABS plastic
ABS basket with divider each drawer
Round rubber bumper
4 inch casters (Two with breaks)



UTILITY CART (RECESSED SHELF)

Model	Size (L×D×H / mm)
KS-UC2S-S	650 × 450 × 900
KS-UC2S-M	750 × 450 × 900
KS-UC2S-L	950 × 450 × 900

Stainless steel 304 construction
3 side guard rail on each shelf
Two recessed shelves
4 inch casters (Two with breaks)
Round rubber bumper
Requested size available

UTILITY CART (RECESSED SHELF)

Model	Size (L×D×H / mm)
KS-UC3S-S	650 × 450 × 900
KS-UC3S-M	750 × 450 × 900
KS-UC3S-L	950 × 450 × 900

Stainless steel 304 construction
3 side guard rail on each shelf
Three recessed shelves
4 inch casters (Two with breaks)
Round rubber bumper
Requested size available

DRESSING CART

Model	Size (L×D×H / mm)
KS-DC2S-S	650 × 450 × 900
KS-DC2S-M	750 × 450 × 900
KS-DC2S-L	950 × 450 × 900

Stainless steel 304 construction
3 side guard rail on each shelf
Two flat shelves
4 inch casters (Two with breaks)
Round rubber bumper
Requested size availableA



UTILITY CART (A.B.S DRAWER)

Model	Size (L×D×H / mm)
KS-AUC2S-H2D	755 × 500 × 950

Stainless steel 304 construction
3 side guard rail on each shelf
Two flat shelves
Drawer front : ABS plastic
Aluminum divider each drawer
(Without removable ABS basket)
Round rubber bumper
4 inch casters (Two with breaks)

UTILITY CART (A.B.S DRAWER)

Model	Size (L×D×H / mm)
KS-AUC3S-H2D	755 × 500 × 950

Stainless steel 304 construction
3 side guard rail on each shelf
Three flat shelves
Drawer front : ABS plastic
Aluminum divider each drawer
(Without removable ABS basket)
Round rubber bumper
4 inch casters (Two with breaks)

UTILITY CART (A.B.S DRAWER)

Model	Size (L×D×H / mm)
KS-AUC2S-H4D	755 × 500 × 950

Stainless steel 304 construction
3 side guard rail on each shelf
Two flat shelves
Drawer front : ABS plastic
Aluminum divider each drawer
(Without removable ABS basket)
Round rubber bumper
4 inch casters (Two with breaks)



DRESSING CART

Model	Size (L×D×H / mm)
KS-DC3S-S	650 × 450 × 900
KS-DC3S-M	750 × 450 × 900
KS-DC3S-L	950 × 450 × 900

Stainless steel 304 construction
3 side guard rail on each shelf
Three flat shelves
4 inch casters (Two with breaks)
Round rubber bumper
Requested size available

DRESSING CART (BENDING TYPE)

Model	Size (L×D×H / mm)
KS-DC2SH-S	650 × 450 × 950
KS-DC2SH-M	750 × 450 × 950
KS-DC2SH-L	950 × 450 × 950

Stainless steel 304 construction
3 side guard rail on each shelf
Handle at both side
Two flat shelves
4 inch casters (Two with breaks)
Round rubber bumper
Requested size available

DRESSING CART (BENDING TYPE)

Model	Size (L×D×H / mm)
KS-DC3SH-S	650 × 450 × 950
KS-DC3SH-M	750 × 450 × 950
KS-DC3SH-L	950 × 450 × 950

Stainless steel 304 construction
3 side guard rail on each shelf
Handle at both side
Three flat shelves
4 inch casters (Two with breaks)
Round rubber bumper
Requested size available



TREATMENT CART

Model	Size (L×D×H / mm)
KS-TC2S1-S	550×450×900
KS-TC2S1-M	650×450×900

Stainless steel 304 construction
3 side guard rail on each shelf
Two flat shelves
One SUS Drawer
4 inch casters (Two with breaks)
Round rubber bumper
Additional shelf or drawer available



TREATMENT CART

Model	Size (L×D×H / mm)
KS-TC2S2-S	550×450×900
KS-TC2S2-M	650×450×900

Stainless steel 304 construction
3 side guard rail on each shelf
Two flat shelves
Two SUS Drawers
4 inch casters (Two with breaks)
Round rubber bumper
Additional shelf or drawer available



TREATMENT CART

Model	Size (L×D×H / mm)
KS-TC2S2-S	750×450×900
KS-TC2S2-M	850×450×900
KS-TC2S2-L	950×450×900

Stainless steel 304 construction
3 side guard rail on each shelf
Two flat shelves
Two SUS Drawers side by side
4 inch casters (Two with breaks)
Round rubber bumper
Additional shelf or drawer available



TREATMENT CART

Model	Size (L×D×H / mm)
KS-TC3S2-S	750×450×900
KS-TC3S2-M	850×450×900
KS-TC3S2-L	950×450×900

Stainless steel 304 construction
3 side guard rail on each shelf
Three flat shelves
Two SUS Drawers side by side
4 inch casters (Two with breaks)
Round rubber bumper
Additional shelf or drawer available



TREATMENT CART (BANDING TYPE)

Model	Size (L×D×H / mm)
KS-TC2S2H-S	750×450×950
KS-TC2S2H-M	850×450×950
KS-TC2S2H-L	950×450×950

Stainless steel 304 construction
3 side guard rail on each shelf
Handle at both side
Two flat shelves
Two SUS Drawers side by side
4 inch casters (Two with breaks)
Round rubber bumper
Additional shelf or drawer available



R-ABS



TREATMENT CART (BANDING TYPE)

Model	Size (L×D×H / mm)
KS-TC3S2H-S	750×450×950
KS-TC3S2H-M	850×450×950
KS-TC3S2H-L	950×450×950

Stainless steel 304 construction
3 side guard rail on each shelf
Handle at both side
Three flat shelves
Two SUS Drawers side by side
4 inch casters (Two with breaks)
Round rubber bumper
Additional shelf or drawer available




SIDE TABLE

Model	Size (L×D×H / mm)	Remark
KS-ST-01	800×350×775	2" Caster, Wood, Steel
KS-ST-02	800×350×775	1" Caster, H45



SIDE TABLE (C-ARM)

Model	Size (L×D×H / mm)	Remark
KS-ST-03	700×400×775	BAKE 12T
KS-ST-05	700×400×775	BAKE 12T



URINE CONTAINER RACK (HOLDER TYPE)

Model	Size (L×D×H / mm)	Remark
KS-UCR01	100×220×260	1 Holder
KS-UCR02	200×220×260	2 Holder
KS-UCR03	300×220×260	3 Holder

All Stainless Steel 304 Construction
Pipe Frame & Joint Type



FECES PEN RACK (HOLDER TYPE)

Model	Size (L×D×H / mm)	Remark
KS-FPR01	100×220×330	1 Holder
KS-FPR02	200×220×330	2 Holder
KS-FPR03	300×220×330	3 Holder

All Stainless Steel 304 Construction
Pipe Frame & Joint Type



SIDE TABLE (ARM SUPPORT)

Model	Size (L×D×H / mm)	Remark
KS-ST-04	460×570×775	mattress 600x150x40



BED PEN RACK

Model	Size (L×D×H / mm)	Remark
KS-BPR03	500×420×800	3 Holder
KS-BPR05	850×420×800	5 Holder
KS-BPR08	1100×420×800	8 Holder
KS-BPR10	1400×420×800	10 Holder

All Stainless Steel 304 Construction
Drain Hole & Tray
Wall Mounted Type



LINEN HAMPER CART

Model	Size (L×D×H / mm)	Remark
KS-LHCF	600×600×850	—
KS-LHCF-1	900×600×850	—
KS-LHCF-2	1200×600×850	—

Stainless Steel Pipe 304 Frame
4" Swivel Caster 4EA
Accessories : Hamper bag 1EA
Option : Water proof hamper bag available
(Color : Orange or Green)



LINEN HAMPER CART

Model	Size (L×D×H / mm)	Remark
KS-LHCR	φ600×850	—

Stainless Steel Pipe 304 Frame
4" Swivel Caster 3EA
Accessories : Hamper bag 1EA



LINEN HAMPER CART (PEDAL TYPE)

Model	Size (L×D×H / mm)	Remark
—	OPTION	—



LINEN HAMPER CART

Model	Size (L×D×H / mm)	Remark
KS-LHCD	1000×500×850	Double

Stainless Steel Pipe 304 Frame
Stainless Steel Shelf Plate
Accessories : Hamper bag 2EA



LINEN HAMPER CART (PEDAL TYPE)

Model	Size (L×D×H / mm)	Remark
KS-LHCF-C	500×500×850	—
KS-LHCF-C6	600×600×850	—

OPTION-Waterproof (orange, green)



LINEN HAMPER BAG (ORANGE/GREEN/IVORY)

Model	Size (L×D×H / mm)	Remark
KS-LHB-01	500×500	—
KS-LHB-02	600×600	—



STEP FOOT STOOL

Model	Size (L×D×H / mm)	Remark
KS-FS1-1	430×300×120	One Step
KS-FS1-2	430×300×180	One Step
KS-FS1-3	430×300×250	One Step

All Stainless Steel 304 Construction
Non Slip rubber top surface
Non Slip rubber Parking



STEP FOOT STOOL

Model	Size (L×D×H / mm)	Remark
KS-FS2-1	430×500×500	Two Step
KS-FS2-2	700×500×500	Two Step

All Stainless Steel 304 Construction
Non Slip rubber top surface
Non Slip rubber Parking



STEP FOOT STOOL (BENDING TYPE)

Model	Size (L×D×H / mm)	Remark
KS-FS-3	430×500×470	Two Step

All Stainless Steel 304 Construction



STEP FOOT STOOL

Model	Size (L×D×H / mm)	Remark
KS-FS3-H	700×500×900	Two Step

All Stainless Steel 304 Construction
Non Slip rubber top surface
Non Slip rubber Parking



CHART HOLDER CART

Model	Size (L×D×H / mm)	Remark
KS-CHC530	530×335×1000	40매
KS-CHC630	630×335×1000	60매
KS-CHC830	830×335×1000	80매

All Stainless Steel 304 Construction
4" Swivel Caster 4EA



CHART HOLDER CART

Model	Size (L×D×H / mm)	Remark
KS-CHC530S	530×335×1200	40매
KS-CHC630S	630×335×1200	60매
KS-CHC830S	830×335×1200	80매

All Stainless Steel 304 Construction
4" Swivel Caster 4EA



BASKET CART

Model	Size (L×D×H / mm)	Remark
KS-BC1	600×450×900	—



BASKET CART

Model	Size (L×D×H / mm)	Remark
KS-BC2	600×450×900	—



OXYGEN TANK CART

Model	Size (L×D×H / mm)	Remark
KS-OTC20	For 20kg	10L

SUS



OXYGEN TANK CART

Model	Size (L×D×H / mm)	Remark
KS-OTC10-A	—	φ120
KS-OTC10-B	—	φ160

All SUS



BASKET CART

Model	Size (L×D×H / mm)	Remark
KS-BC4	φ120×440	5.5L



BASKET CART

Model	Size (L×D×H / mm)	Remark
KS-BC3	700×500×950	—



OXYGEN HOLDER

Model	Size (L×D×H / mm)	Remark
KS-OTC6-115	390X260X230	—



OXYGEN TANK CART

Model	Size (L×D×H / mm)	Remark
KS-OTC40	For 40kg	20L

Steel



OXYGEN TANK HOLDER

Model	Size (L×D×H / mm)	Remark
KS-OTH120	φ120×440	5.5L



OXYGEN TANK HOLDER


Model	Size (L×D×H / mm)	Remark
KS-OTH160	φ160×370	5L



OXYGEN TANK CART

Model	Size (L×D×H / mm)	Remark
KS-OTC6-116	490X315X900	—

SUS



SCREEN

Model	Size (L×D×H / mm)	Remark
KS-1PS	900×440×1800	1 Panel

Stainless steel 304 pipe frame
Easy to change fabric screen
2 inch swivel casters



SCREEN

Model	Size (L×D×H / mm)	Remark
KS-2PS	1800×400×1800	2 Panel

Stainless steel 304 pipe frame
Easy to change fabric screen
2 inch swivel casters
Foldable each screen side



SCREEN

Model	Size (L×D×H / mm)	Remark
KS-3PS	2700×400×1800	3 Panel

Stainless steel 304 pipe frame
Easy to change fabric screen
2 inch swivel casters
Foldable each screen side



UTILITY CART

Model	Size (L×D×H / mm)	Remark
KS-UC3S-W3	1030X450X870	



EMERGENCY CART (K TYPE)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-KB	740X575X985(1055)	BLUE

Top shelf : ABS plastic with SUS 304
Side & Back panel : Aluminum
Guard rail : Steel (Color : Grey)
Column : Aluminum (Color : Grey)
Drawer front : Powder Coated Steel
Drawer with removable ABS basket
Drawer rail : Auto return available upon request
Total key locking
Caster : 5 inch (Two with brake)



EMERGENCY CART (K TYPE)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-KY	740X575X985(1055)	YELLOW

Top shelf : ABS plastic with SUS 304
Side & Back panel : Aluminum
Guard rail : Steel (Color : Grey)
Column : Aluminum (Color : Grey)
Drawer front : Powder Coated Steel
Drawer with removable ABS basket
Drawer rail : Auto return available upon request
Total key locking
Caster : 5 inch (Two with brake)



EMERGENCY CART (K TYPE)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-K7	740X575X985(1055)	5 COLORS

Top shelf : ABS plastic with SUS 304
Side & Back panel : Aluminum
Guard rail : Steel (Color : Grey)
Column : Aluminum (Color : Grey)
Drawer front : Powder Coated Steel
Drawer with removable ABS basket
Drawer rail : Auto return available upon request
Total key locking
Caster : 5 inch (Two with brake)



EMERGENCY CART (K TYPE)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-KSTS	740X575X985(1055)	SUS SHELF

Top shelf : Stainless Steel 304
Side & Back panel : Aluminum
Guard rail : SUS304
Column : Aluminum (Color : Grey)
Drawer front : ABS plastic
Drawer with removable ABS basket
Drawer rail : Auto return available upon request
Total key locking
Caster : 5 inch (Two with brake)



EMERGENCY CART (K TYPE)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-STY	660X510X1000(1070)	YELLOW

Top shelf : Stainless Steel 304
Body frame : Powder Coated Steel
Guard rail : SUS304
Drawer front : Powder Coated Steel
Drawer with removable ABS basket
Caster : 5 inch (Two with brake)



EMERGENCY CART (K TYPE)

Model	Size (L×D×H / mm)	Remark
KS-EMC5D-STP	660X510X1000(1070)	PINK

Top shelf : Stainless Steel 304
Body frame : Powder Coated Steel
Guard rail : SUS304
Drawer front : Powder Coated Steel
Drawer with removable ABS basket
Caster : 5 inch (Two with brake)