


350 SERIES




ENDO 350 CART

Model	Size (L×D×H / mm)	Remark
KS-ES350-480S	640×670×1220	slide type keyboard

모니터암 1ea / 내시경 행거 / 220×6구

OPTION



▪ 키보드 슬라이드



▪ 모니터 암(S암)
KS-SARM



▪ 내시경 행거(ㄱ타입)



▪ 사이드 홀더
KS-PIN-BL



▪ 센터 홀더



▪ 카메라 보호대
KS-SH100



ENDO 350 (PREMIUM)

Model	Size (L×D×H / mm)	Remark
KS-ES350-P	640×670×1220	-

모니터암 2ea / 내시경 행거 / 220×10구



ENDOSCOPE CART

Model	Size (L×D×H / mm)	Remark
KS-CS123	525×470	19 봉

VESA-75x75 / 100x100 / 9kg~19kg



ASSIST MONITOR CART

Model	Size (L×D×H / mm)	Remark
KS-HDMC-01	φ 720×400~1900	27" 이하


기본 2kg~8kg / 고하중 8kg~13kg




ASSIST MONITOR CART

Model	Size (L×D×H / mm)	Remark
KS-HDMC-N	710X700X1750~2130	40~65 inch

450 SERIES



Detail image




ENDO 450 CART

Model	Size (L×D×H / mm)	Remark
KS-ES450-480S-OL	630×660×1180	—

Shelf: 480X550 Weight: 52kg

OPTION



키보드 슬라이드



모니터 암 (SARM)
KS-SARM




내시경 행거(그타입)



사이드 홀더
KS-PIN-BL



센터 홀더



카메라 보호대
KS-SH100



ARTHROSCOPE CART

Model	Size (L×D×H / mm)	Remark
KS-AS350-480S	630×670×1300	—

Shelf: 480X545

SPEC: 5" Caster, Lcd Arm(s), 220V x 6, Keyboard Slide

ADD: IV Hnager, Oxygen tank Holder



ARTHROSCOPE CART

Model	Size (L×D×H / mm)	Remark
KS-AS350-480S	630×670×1300	—

Shelf: 480X545

SPEC: 5" Caster, Lcd Arm(s), 220V x 6, Keyboard Slide

ADD: IV Hnager, Oxygen tank Holder



ENDOSCOPE CART

Model	Size (L×D×H / mm)	Remark
KS-ES600-480G	650×795×1450	G.DOOR



ARTHROSCOPE CART

Model	Size (L×D×H / mm)	Remark
KS-AS450-480-4	630×650×1255	4 Plate
KS-AS450-480-5	630×650×1355	5 Plate

Shelf: 480X530

SPEC: 4" Caster, Lcd Arm(s), 220V x 6

ADD: IV Hnager, Keyboard Slide, Oxygen tank Holder

Operating Room
Appliances

Nursing
Appliances

Consulting Room
Appliances

Supply Room
Appliances

Hospital Bed

Carts

Monitor Carts

IV Pole Stand

Monitor Arm

Wall Mount


Pole Mount

Accessory

W350 SERIES



Without pump hook



LAPAROSCOPE CART LAPA 350 CART

Model	Size (L×D×H / mm)	Remark
KS-LS350-480	640×670×1300	—

Shelf: 480X545
SPEC: 5" Caster, Oxygen tank Holder, Drawer, Lcd Arm(s), Pump Hanger
ADD: Keyboard Slide, IV Hanger

OPTION



키보드 슬라이드(중판)



모니터암(S암)
KS-SARM




산소통받침대(W180)



사이드 홀더
KS-PIN-BL



센터 홀더



카메라 보호대
KS-SH100



펌프걸이



LAPAROSCOPE CART

Model	Size (L×D×H / mm)	Remark
KS-LS350-480S-LOW	640×670×1300	—

Shelf: 480X545
SPEC: 5" Caster, Oxygen tank Holder, Drawer, Lcd Arm(s), Keyboard Slide, Pump Hanger, IV Hanger, 220V x 6



ARTHROSCOPE CART

Model	Size (L×D×H / mm)	Remark
KS-AS-LVT-A	520×700×1345	—

Shelf: 370X565
SPEC: 5" Caster, Stand Arm
ADD: Joint Arm, Extension Arm, Keyboard Slide



LAPAROSCOPE CART

Model	Size (L×D×H / mm)	Remark
KS-LS450-480-BB	630×670×1300	—

Shelf: 480X545
4 Caster, Oxygen Tank Holder (2EA) Back



ARTHROSCOPE CART

Model	Size (L×D×H / mm)	Remark
KS-AS-AR	550×600×1345	—

Shelf: 400X440
SPEC: 5" Caster, Keyboard Slide, Joint Arm, Extension Arm, Wire Basker



ARTHROSCOPE CART

Model	Size (L×D×H / mm)	Remark
KS-AS-LVT-B	520×700×1345	—

Shelf: 370X565
SPEC: 5" Caster, Stand Arm, Keyboard Slide
ADD: Joint Arm, Extension Arm

Operating Room
Appliances

Nursing
Appliances

Consulting Room
Appliances

Supply Room
Appliances

Hospital Bed

Carts

Monitor Carts

IV Pole Stand


Monitor Arm

Wall Mount

Pole Mount

Accessory


W350 SERIES



E.N.T CART

Model	Size (L×D×H / mm)	Remark
KS-ENT350-480	630×670×1090	—


Shelf: 480X445
SPEC: 5" Caster, Lcd Arm(s), Acrylic Hanger, 220V×6
ADD: Keyboard Slide, Drawer



E.N.T CART

Model	Size (L×D×H / mm)	Remark
KS-ENT450-480	630×560×1045	—

Shelf: 480X445
SPEC: 4" Caster, Acrylic Hanger, 220V×6
ADD: Keyboard Slide, Drawer, Lcd Arm(s).



PC CART

Model	Size (L×D×H / mm)	Remark
KS-PC350	680×620×835	—

Shelf: 520X495
SPEC: 5" Caster, Lcd Arm(s), Keyboard Slide, 220V x 6
ADD: Joint Arm, Drawer




PC CART

Model	Size (L×D×H / mm)	Remark
KS-PC350D	680×620×832	—


Shelf: 520X495
SPEC: 5" Caster, Keyboard Slide, 220V x 6
ADD: Joint Arm, Lcd Arm(s).



소독관4구 거치대
KS-SHB-4
KS-SHB-2
KS-SHB-1



ENDO KEEP



소독관스탠드 카트(1구)
KS-RS-SHB1



E.N.T CART

Model	Size (L×D×H / mm)	Remark
KS-ENT450-480W	630×560×1045	—

Shelf: 480X440
SPEC: 4" Caster, 220V x 6
ADD: Lcd Arm(s), Keyboard Slide, Drawer



EMR CART

Model	Size (L×D×H / mm)	Remark
KS-EMR0355	480X580X850	—




PC CART

Model	Size (L×D×H / mm)	Remark
KS-PC350-22D	—	—



PC CART

Model	Size (L×D×H / mm)	Remark
KS-PC350-22D3	—	—



PC CART

Model	Size (L×D×H / mm)	Remark
KS-PC450	670×620×800	—

Shelf : 520X495
SPEC : 4" Caster, Lcd Arm(s), Keyboard Slide, 220V×6
ADD : Joint Arm, Drawer



PC CART

Model	Size (L×D×H / mm)	Remark
KS-PC350S	680×620×1090	—

Shelf : 520X495
SPEC : 5" Caster, Keyboard Slide, 220V×6
ADD : Lcd Arm(s), Drawer Plate

W350 SERIES



PC CART

Model	Size (L×D×H / mm)	Remark
KS-PC450-SDM	670×600×1010	—

Shelf: 520X415
SPEC : 3" Caster, Lcd Arm(s), Mouse Plate, Drawer 220V x 4



PC CART

Model	Size (L×D×H / mm)	Remark
KS-PCSU-D1	500×475×1045	—

4" Caster, Lcd Arm(s), Drawer SUS 304



PC CART

Model	Size (L×D×H / mm)	Remark
KS-PC350-DM	620×570×1070	—

Shelf : 480X460
SPEC : 4" Caster, Mouse Plate, PC Bracket, 220V×6
ADD : Lcd Arm(s), Joint Arm



PC CART

Model	Size (L×D×H / mm)	Remark
KS-PC450-DM	610×520×960	—

Shelf : 470X450
SPEC : 4" Caster, Mouse Plate, PC Bracket, 220V×6
ADD : Lcd Arm(s), Joint Arm



PC CART

Model	Size (L×D×H / mm)	Remark
KS-PCSU-DS	500×350×1063	—

4" Caster, Lcd Arm(s), Keyboard Slide, Drawer, SUS 304



PC CART

Model	Size (L×D×H / mm)	Remark
KS-PC-KC01	510×575×885	—

Shelf : 490X460
3" Caster, Keyboard Slide



BOBBY CART

Model	Size (L×D×H / mm)	Remark
KS-FT10	470×650×1000	—



EMR CART

Model	Size (L×D×H / mm)	Remark
제품명	L×D×H	—



CCTV CART

Model	Size (L×D×H / mm)	Remark
제품명	L×D×H	—



제품명

Model	Size (L×D×H / mm)	Remark
제품명	L×D×H	—



NOTEBOOK CART

Model	Size (L×D×H / mm)	Remark
KS-PC-YW01	600×575×880	—

Shelf: 600X515
4" Caster, Drawer, Handle, 220V×6



EEG CART

Model	Size (L×D×H / mm)	Remark
KS-PC-KW01	560×580×830	—

Shelf: 500X500
5" Caster, Handle, Keyboard Slide, 220V×6



MULTI CART

Model	Size (L×D×H / mm)	Remark
KS-MT350-D2	590×520×815	—

Shelf: 440X460
SPEC : 4" Caster, 4 Handle, Drawer
ADD : Keyboard Slide



MULTI CART

Model	Size (L×D×H / mm)	Remark
KS-MT500-L1	640×520×860	—

Shelf: 500X463
SPEC : 5" Caster, 4 Handle
ADD : Keyboard Slide, Drawer, Basket, Plate



MULTI CART

Model	Size (L×D×H / mm)	Remark
KS-MT350	590×520×815	—

Shelf: 440X460
SPEC : 5" Caster, 4Handle, Keyboard Slide, Basket, Drawer, Plate



MULTI CART

Model	Size (L×D×H / mm)	Remark
KS-MT350-B	590×520×815	—

Shelf: 440X460
SPEC : 5" Caster, 4 Handle, Basket
ADD : Keyboard Slide, Drawer



MULTI CART

Model	Size (L×D×H / mm)	Remark
KS-MT500-P1	640×520×860	—

Shelf: 500X463
SPEC : 5" Caster, 4 Handle
ADD : Keyboard Slide, Drawer, Basket



MULTI CART


Model	Size (L×D×H / mm)	Remark
KS-MT450-D3	670×620×800	—

Shelf: 520X495
SPEC : 4" Caster, Drawer, 220V x 6
ADD : Lcd Arm(s), Joint Arm



PC CART(DOUBLE)

Model	Size (L×D×H / mm)	Remark
KS-PC350-TB2	—	—



PC CART

Model	Size (L×D×H / mm)	Remark
KS-EMR0354-S	510×605×1090	—
KS-EMR0354-L	510×605×822	—



MULTI CART

Model	Size (L×D×H / mm)	Remark
KS-MT4-3	—	—



MULTI CART

Model	Size (L×D×H / mm)	Remark
KS-MT4-4	—	—



MULTI CART

Model	Size (L×D×H / mm)	Remark
KS-MT450-S-B	580×470×730	—

Shelf : 480X430
SPEC : 3" Caster, Keyboard Slide
ADD : Drawer, Plate



MULTI CART

Model	Size (L×D×H / mm)	Remark
KS-MT350-US	680×620×835	—

Shelf : 520X610
SPEC : 5" Caster, Drawer, Gel Bracket
ADD : Lcd Arm(s), Joint Arm



BOVIE CART

Model	Size (L×D×H / mm)	Remark
KS-BV-BW	610×570×1000	—

Shelf : 450X300, 450
5" Caster, Drawer, Basket



BOVIE CART

Model	Size (L×D×H / mm)	Remark
KS-BV-VL-UC	500×560×880	—

Shelf : 470X495
5" Caster, The Upper Shelf, Oxygen Tank Holder



MULTI CART

Model	Size (L×D×H / mm)	Remark
KS-MT500-DB	640×460×800	—

Shelf : 500X400
3" Caster, Drawer, Bracket, Basket



BOVIE CART

Model	Size (L×D×H / mm)	Remark
KS-BV450-D	540×520×840	—

Shelf : 400X500
SPEC : 4" Caster, Drawer
ADD : Lcd Arm(s), Drawer Basket



BOVIE CART

Model	Size (L×D×H / mm)	Remark
KS-BV-UES40	380×470×815	—

Shelf : 370X420
4" Caster, Steel, Corner Safe Rubber



BLOOD PRESSURE SYSTEM

Model	Size (L×D×H / mm)	Remark
KS-SM-01	460×400×800	9020

SPEC : Adjust Foot Rubber
ADD : Wire Basket, Caster



MULTI CART(CTX)

Model	Size (L×D×H / mm)	Remark
KS-MT4-1	490×515×745	—



MULTI CART(오스테오)

Model	Size (L×D×H / mm)	Remark
KS-MT4-2	412×525×890	—



혈압계 CART

Model	Size (L×D×H / mm)	Remark
KS-SM-9030	460×400×780	—



제품명

Model	Size (L×D×H / mm)	Remark
KS-D50		



EKG CART

Model	Cart Size (L×D×H/mm)	Inner Size (L×D×H/mm)
KS-EG-S	405×325×740	370×285×15
KS-EG-M	530×325×740	490×280×15
KS-EG-L	485×405×740	440×365×15
KS-EG-B	540×465×740	490×420×15

3" Caster, General Hanger, Steel



EKG CART (STANDERD HANGER)

Model	Cart Size (L×D×H/mm)	Inner Size (L×D×H/mm)
KS-EG-S-L	405×325×740	370×285×15
KS-EG-M-L	530×325×740	490×280×15
KS-EG-L-L	485×405×740	440×365×15
KS-EG-B-L	540×465×740	490×420×15

3" Caster, Standard Hanger, Steel



STANDARD HANGER

Model	Size (L×D×H / mm)	Remark
KS-EGH-L	H=610	φ 14 PIPE



JOINT HANGER

Model	Size (L×D×H / mm)	Remark
KS-NEGH-H	H=860	φ 14 PIPE



EKG CART (JOINT HANGER)

Model	Cart Size (L×D×H/mm)	Inner Size (L×D×H/mm)
KS-EG-S-H	405×325×740	370×285×15
KS-EG-M-H	530×325×740	490×280×15
KS-EG-L-H	485×405×740	440×365×15
KS-EG-B-H	540×465×740	490×420×15

3" Caster, Joint Hanger, Steel



EKG CART

Model	Size (L×D×H / mm)	Remark
KS-EG-7402	485×405×820	Hole (x)

SPEC : 3" Caster, Steel
ADD : Wire Basket
Plate : Outside 480x405x25 / Inside 440x365x15



EKG CART

Model	Size (L×D×H / mm)	Remark
KS-SM-02	485×405×820	—

3" Caster, Wire Basket
Plate : Outside 480x405x25 / Inside 440x365x15



EKG CART

Model	Size (L×D×H / mm)	Remark
KS-EG-7402-H	485×405×820	Hole (o)

SPEC : 3" Caster, Steel, Paper supply space
ADD : Wire Basket
Plate : Outside 480x405x25 / Inside 440x365x15



EKG CART

Model	Size (L×D×H / mm)	Remark
KS-EG-1350	300×380×900	—

4" Caster



EKG CART

Model	Size (L×D×H / mm)	Remark
KS-EG-DP	320×350×878	—

3" Caster, Keyboard Plate, Basket



NEBULIZER

Model	Size (L×D×H / mm)	Remark
KS-NE-01	415×450×875	—

SUS 304
Nebulizer



BEAUTY

Model	Size (L×D×H / mm)	Remark
KS-EG450	310×370×733	—

3" Caster



MEB

Model	Size (L×D×H / mm)	Remark
KS-EG-MEB	500×608×820	—

4" Caster, Mouse Plate, Lcd Arm



BEAUTY-B850



Model	Size (L×D×H / mm)	Remark
KS-GEG000I	600×1370	—

3" Caster



FOCUS CART

Model	Size (L×D×H / mm)	Remark
KS-PKS100	660×650×H875	—



ASSIST MONITOR CART (MOTOR)


Model	Size (L×D×H / mm)	Remark
KS-HDMC-02	φ 630×1430~1930	—

2,5" Caster, Steel, Aluminium, Plastic, Motor
Using 220V

Easily move the monitor up and down.
You can adjust the height using the remote controll
You can adjust the VESA HOLE depending on the type of monitor.
220 volts of electricity



Remote control



VESA hole



Adjust hole spacing



Handle



Adjust hole spacing